

ForestGEO - sampling protocol for wood density

Norbert Kunert, May 2019

1. General information

Per species, at least three tree individuals should be sampled. You can use the botanical material collected for the leaf trait analysis.

2.1. List of material needed

- Pruner, or other tools to cut branches
- Small receptacle with low weight (total weight including water should be within the range of your scale)
- Flagging tape to mark branches
- Waterproof permanent marker
- Razor blades or knife
- Dissection needle
- Scale (precision 0.01g)
- Water
- Paper towels
- Oven to dry samples


2.2. Preparation of sample

- Cut a segment of 2-3cm from each sampled tree (you can use the part you re-cut under water from the leaf trait measurements).
- Mark segments with flagging tape (tree species and individual).
- Put in water over night.


3. Wood density measurement

3.1 Fresh sample

- Take sample out of the water and dry with paper towel.
- Measure the stem segment diameter in the middle.
- Measure the length of the segment. Weigh the stem segment.
- Remove the bark with knife or razor blades.
- Measure the stem diameter without bark
- Split the stem segment along the pith and remove all soft tissue.
- Weigh the two halves of the segments again (you can weigh both halves together).



- Place receptacle with water on the scale and press Tara.
- Pin the two halves of the segment on the dissection needle and submerge in the water. Note: submerge the halves one by one and add up the two measurements.
- Dry the samples in an oven at 60°C for at least 48 hours.



3.2 Dry sample

- Weigh the dried wood again.
- Wood density is calculated as dry mass divided by volume of displaced water (g m^{-3})